



**Mobile functional furniture
for the healthcare sector**



Catalogue valid from 2020

novocal – care in motion.

As a manufacturer of functional mobile furniture for Health Care, it is our desire to optimize and simplify the working life of nursing staff in hospitals, nursing care and retirement homes. Naturally, your patients are at the centre of your business, and our products allow staff to move from location to location quickly, smoothly and efficiently maximizing your care potential and reducing stress.

We are a technically broad and experienced team based in Ostrhauderfehn in the North of Germany. As a customer you will receive all your bespoke service needs from a single source: advice and quoting, processing and production as well as delivery and instruction.

All our products are made according to individual requirements with service, quality and fine detail our top priority.

Our product range includes functional furniture in the field of laundry logistics and waste management as well as general and special care trolleys, e.g. MRSA, endoscopy or module trolleys. In addition, we also offer all kinds of transport equipment, storage systems and components for sterile surgical areas. Numerous accessories are available to customize your products to suit your individual needs.

Within this catalogue you will find examples of both our wide ranging and individual assortment of products. If you have any special requests or enquiries, please do not hesitate to contact us. We will happily find a solution to any of your requirements. Our customers trust in us to deliver our promises.

Together we can set care in motion.



*Our promise:
Everything from a single source.
Reliable, sustainable, individual.
Your Novocal team*

Laundry Logistics

	Article Number	Catalogue Page
Linen Collectors Stainless Steel		
Wall mounting collector	NW10	11
„Compact“ 1-fold / „Compact plus“ 1-fold with push bar	NW71/NW81	11
„Compact“ 2-fold / „Compact plus“ 2-fold with push bar	NW72/NW82	12
„Compact“ 3-fold	NW73	12
„Compact plus“ 3-fold with push bar	NW83	13
„Compact“ 4-fold / „Compact plus“ 4-fold with push bar	NW74/NW84	13
„Compact“ 5-fold / „Compact plus“ 5-fold with push bar	NW75/NW85	13
„Compact“ 5-fold with centralised lid operating system	NW25-G	13
„Compact plus“ 3-fold to 5-fold with push bar	NW43F/NW44/NW45	14
Combination collector 2-fold / 3-fold with push bar	NW42-KS70/NW43-KS70	14
Combination collector 1-fold to 3-fold with push bar	NW81-KS70/NW82-KS70/NW83-KS70	15
Combination collector 4-fold with push bar	NW84-KS70	15
Accessories Linen Collectors		
Plastic lid with hook holder	KDH	16
Plastic lid w. plastic hinge & soft close system (coloured/transp.)	KDS/KDS-TR	16
Pedal device for hinged lid	TVD	16
Adhesive ring	SR400	16
Synthetic base tub for collectors 1-fold to 4-fold	KBW	16
Linen Bags		
Linen bags, different versions	TW1/TW2	17
Linen bags, waterproof and water-soluble	WTW2/HWW	17
Linen Nets		
Linen nets, different versions	N/NF	18
Toggle/Toggle made of rubber	DIBE/GK	18
Stainless Steel Undercarriages		
Stainless steel undercarriage with plastic tub/with push bar	WW1/WW2	19
Laundry Trolleys		
Overview Laundry trolleys	SW500-SW2000	20-21
Possible arrangements of adhesive rings		22
Laundry trolleys of the series SW500	SW500	23
Laundry trolleys of the series SW600	SW600	24
Laundry trolleys of the series SW700	SW700	25
Laundry trolleys of the series SW800	SW800	26
Laundry trolleys of the series SW900	SW900	27
Laundry trolleys of the series SW1000	SW1000	28
Laundry trolleys of the series SW2000	SW2000	29

Patient Care and Nursing Trolleys

PAW		
Overview Patient Care and Nursing Trolleys	PAW3000-PAW7000	30-31
Individual production		32
Patient care and nursing trolleys of the series PAW3000	PAW3000	33
Patient care and nursing trolleys of the series PAW4000	PAW4000	34
Patient care and nursing trolleys of the series PAW5000	PAW5000	35
Patient care and nursing trolleys of the series PAW6000	PAW6000	36
Patient care and nursing trolleys of the series PAW7000	PAW7000	37
Patient care and nursing trolleys of the series PAW8000	PAW8000	38
Patient care and nursing trolleys of the series PAW9000	PAW9000	39

Contents



	Article Number	Catalogue Page
Medicine Trolleys		
Medicine trolleys, various versions	MEDI5001/MEDI6002	40
Accessories		41
Module Trolleys		
Module trolleys series 6000	MOD	42
Module trolleys, open versions	OMW	43
Roll-Up utility trolley	RUWI	43
Modules		
ISO modules ABS, divisible, different sizes	ABS	44
ISO modules PC, divisible, different sizes	PC	45
Wire-netted baskets, different sizes	GK	45
Baby Changing Trolleys		
Baby changing trolleys, various versions	WW	46
System Trolleys		
System trolley "Modula"	NSW/NSWK	47
Endoscopy Trolley		
Endoscopy trolley	ENDO61000	48
Document Trolleys		
File and ward round trolley	NKV	49
Ergonomic documentation trolley, various versions	EDW	50
Catheter Trolleys		
Catheter trolley, various versions	KW1/KW2	51
Functional Trolleys with synthetic Surfaces		
Phoenixcare Cetos, various versions	CETOS	52-53
Patient trolley	PATW	54
Computer and IT Trolleys		
Computer Trolleys		
Computer trolley I, various versions	CW4000	55
Computer trolley II, various versions	CW5000	56
Computer trolley III, various versions	CW6000	57
IT trolley I	EDVW5400	58
IT trolley II	EDVW6102	59
Laptop Trolleys		
Laptop trolley (lowerable laptop compartment), various versions	LAW4000A	60
Laptop trolley (sideward-extendable laptop compartment), various versions	LAW4000S	60
Laptop trolley, height adjustable	LAW4300H	61
iPad Trolley		
iPad trolley	IPW4400	62
Power Supply Systems		
Power Supply System	PSS	63

MRSA / ISO Trolleys

	Article Number	Catalogue Page
Disinfection Stations		
DESI, mobile/floor model	DESI1/DESI2	64
Desicenter, stainless steel	DESICENTER-VA	65
Desicenter, coloured	DESICENTER-F	65
UV-C disinfection unit	UEKG01	66-67
MRSA trolley		
MRSA trolley with closed back and racks	MRSA1F	68
MRSA trolley with closed back and girder support system	MRSA1M	68
MRSA trolley with stainless steel back and racks	MRSA2F	69
MRSA trolley with stainless steel back and girder support system	MRSA2M	69

Accessories

Sideways fold-away working surface	SPA1	71
Extendable writing surface	SPA2	71
Sideways height-adjustable working surface (installation on the left/right)	SPA1H-L/SPA1H-R	71
Stainless steel plating, different sizes	EAK	71
ABS synthetic top, different sizes	ABSTOP	71
Tubular frame, different sizes	RGE	71
Side compartments	SF50	72
Girder support systems for ISO modules	MTW	72
Dividers for drawers - lengthwise/transverse	LTS/QTS	72
Surgical instrument container, incl. stainless steel fitting	IAB	72
Swab dispenser, euro-rack model	TSE	72
Medibox (septobox / hallow needle box), 1,5 l	MB	72
Medibox mounting, various versions	MBH	72
Medibox mounting (set)	MBHS	72
Sliding handle, for installation on the left/right	BSL/BSR	72
Disinfectant dispenser, different sizes and versions	DMS	73
Disinfectant dispenser, 500 ml, automatic	DMSA	73
Drip pan, different sizes and versions	TS	73
Disinfectant closure plate, lockable	DMSV	73
Stainless steel pump with dosing limitation for disinfectant dispenser	EPD	73
Euro-rack, stainless steel, single-row/double-row	ES1/ES2	74
Euro-rack, stainless steel, height adjustable, single-row/double-row	ES1H/ES2H	74
Euro-rack, stainless steel, for lateral mounting	ESS	74
T-mounting	TH	74
F-mounting, for installation on the left/right	FHL/FHR	74
Stainless steel container for Euro-rack	EBHES	74
Implement basket	UTBES	74
Glove box mounting, various versions	HSB	75
Bins, various versions and sizes	SCHUETTEN	76
Wire-netted baskets	GK	76
Coat hooks, stainless steel	PHE	76
Coat hooks, plastic	PHK	76
Stainless steel kidney dish	NRS01	76
Catheter tubes, different sizes and versions	KK	76
Mounting for mouth protection box	MUBH	76
Automatic sensor waste bin, 8 l	KB8A	77
Automatic sensor waste bin, 8 l, incl. stainless steel fitting	KB8A-E	77
Synthetic waste bin, 9 l for PAW and SW	KB9	77
Synthetic waste bin, 9 l, incl. stainless steel fitting	KB9-E	77
Synthetic waste bin, 15 l for PAW and SW	KB15	77
Synthetic waste bin, 15 l, incl. stainless steel fitting	KB15-E	77
Hygienic waste bin SLIDE, various sizes and versions	SLIDE	77

	Article Number	Catalogue Page
Cylindrical twist lock	ZY	78
Combination lock in two versions	ZS	78
Electronical numeric code pad	BSML	78
Transponder lock	TPS01	78
Monitor mounting	MH	78
iPad mounting	IPH	78
Furniture safe, various versions	MT1/MT2	78
Drip stand holder / IV-pole holder	IFH	79
Oxygen-cylinder holder	SFH01	79
Master adjusting device for patient care and nursing trolleys	ZFS	79
Magnetic holder	HMF	79
Door holder, central	DFST	79
Wheels, various versions		80
Shelf for defibrillator	ADEF	80
Reanimation panel	RP01	80
help-System/fire extinguisher system		80

Waste Collectors

Waste Collector, Stainless Steel

Waste collector „Compact“ 1-4-fold, stainless steel	NA	81
---	----	----

Waste Collector „Genius“

Waste collector „Genius“, floor model, 40 l/70 l	S40/S70	82
Waste collector „Genius“, mobile, 40 l/70 l	S40L/S70L	82
Waste collector „Genius“, mobile, 70 l, 2-fold	S70L-2	82
Waste collector „Genius“ with braking function	S70L-B	83
Accessories for waste collector „Genius“		83

Waste Collector with Bag Cassette/with Vacuum Unit

Bag cassette, various versions	LGP	84
Collector for bag cassette	NAL	84
Waste sealing unit with closed undercarriage, various versions		85
Waste sealing unit with open undercarriage, various versions		85
Waste collector „Genius Plus“	S70LPLUS	86
Valve bags	VSTP01	86
Vacuum device NOVOVAC	NOVOVAC01	87
Vacuum bags, white/transparent	VAKWS01/VAKTP01	87

Recyclable collector

Recyclable collector for synthetic containers, 50 l, 1-4-fold	GSE	88
Recyclable collectors of the series 2000, various versions	NWE	89
Room recyclable collector, 2-4-fold	ZSE	90
Recyclable collector with coloured synthetic boxes, various versions	WKK	91
Recyclable cabinet of the series 3000, 1-4-fold	ZWK	92
Recyclable cabinet, stainless steel, 1-4-fold	EWK	93
Recyclable cabinet, aluminium, 1-4-fold	NWKA	94

Waste collector with pedal

Pedal dust bin, stainless steel	EE	95
Pedal dust bin, galvanized, various sizes	EL	95
Pedal dust bin, stainless chromium steel, various sizes	EC	95
Closable dustbin, stainless steel, various sizes	GAE	95

Transport and Storage systems

	Article Number	Catalogue Page
Drug Cabinets		
Drug cabinets	MS	96-97
Storage System		
Series FLEXISTORE	FLEXISTORE	98-99
Multi-purpose Trolleys		
Multi-purpose trolley, stainless steel, solid core shelves	SEV	100
Multi-purpose trolley, stainless steel, melamine resin coated shelves	SEM	100
Multi-purpose trolley „Silverline“	S	101
Multi-purpose trolley, synthetic	4091/4093	101
Multi-purpose trolley, stainless steel, closed	SEG	102
Multi-purpose trolley, stainless steel, open	SE	102
Serving Trolleys		
Tea cart	FS	103
Tray Trolley		
Tray trolley, various versions	TW	103
Tray trolley, stainless steel, various versions and sizes	TBW	104
Sorting Trolleys		
Combined sorting trolleys, various versions and sizes	KSR	105
Coat stand, various versions and sizes	KI	105
Rack trolleys, various versions and sizes	RA	106
Sorting Rack Trolley		
Sorting rack trolley, various versions and sizes	SR	107
Standing racks, various versions and sizes	ST	108
Compartment dividers	F	109
Corner posts	P	109
Grid shelves	GB	109
Border grids	RG	109
Textile protective cover	H	110
Pull-out clothing rail	AK	110
Label holder	EK	110
Border poles	RSB	110
Sliding handle for rack trolleys	BSR	111
S-shaped hooks	SH	111
Brackets for border poles	KKL	111
Buffer plates	SSR	111
Directional locking	RIFIX	111
Base plate	FP	111
Transport Trolleys		
Skeleton container trolley	GTR	112
Transport trolley, stainless steel, foldable shelves	TRE91	112
Transport trolley, stainless steel, with / without back, various sizes	TRE	113
Front wall trolley	STE	113
Cabinet Trolley		
Module trolley, closed, various versions	GM	114
Cabinet trolley, aluminium, various versions and sizes	GS	115-117
Spring-base trolley, various sizes	GFW	118

Contents

	Article Number	Catalogue Page
Transport containers, aluminium, water-tight	GTC1	118
Transport container, aluminium, with side panel cut-outs, various sizes	GTC2	118
Transport container, aluminium, hinged lid	GTC350	119
Transport container, aluminium, open	GTW	119
Transport container, aluminium, with foldable front panel	GTWK	119
Transport boxes, aluminium, various sizes	TK	120
Transport container, aluminium, perforated	GTW/L	121

OP Equipment

Transport Trolleys for Sterile Materials

Transport trolleys for sterile materials for inner module frames	CI-STGW	122
Inner frame for transport trolleys	CI-IRA	123
Transport trolleys for sterile materials with bearing rails for modules	CI-STGW	124-125
Transport trolleys for sterile materials with fixed/movable shelves	CI-STGW-FB/-A	126
Rack trolley for sterile materials	CI-STGR	127
Sterile cabinet	CI-STG-SCHRANK	127

OP Equipment

Drip stand/IV-pole	IS	128-129
Accessories for drip stand/IV-pole		130-131
OP stool, various versions	CI-OPH	132
OP tables, various versions	CI-OPB	133
Table for OP instruments, various versions	IE01/CI-IES	134
Working and storage table, mobile and floor model	ATE/CI-ATE	135

Stainless Steel Fabrications

Operation room step	OPT	136
Bowl retainer	SST01	136
Trolley for abdominal bandages	BTW01	136
Hook trolley	HW01	136
Room cleaning trolley	RW01	137
Lab coat stand, stainless steel, various sizes	KI	137
Mangle rack, stainless steel, various sizes	MB	137
Stainless steel mounting for urine bottles, various versions	UF	138
Stainless steel mounting for toilet pails, various versions	TE	138
Stainless steel mounting for bedpans, various versions	SB	138
Stainless steel wall shelf	WSU	138
Trolley for bedpans and urine bottles	TW-US8	138
Shoe rack, various versions	WSR/SRW	139

Seating

ModernArt	140
Relax'nFit	141-142

Index

Index of abbreviations and article numbers in alphabetic sequence	144-145
---	---------

Overview of available colours

					
White (WS)	Vanilla (VL)	Pastel Yellow (PG)	Lemon Yellow (ZG)	Orange (OR)	Orchid (OC)
					
Red (RT)	Lime Green (LI)	Aqua (AQ)	Light Blue (HB)	Steel-blue (SB)	Silver (SI)
					
Light Grey (LG)	Aubergine (AU)	Anthracite (AZ)	Cashmere (CS)	Basalt (BS)	Lava Grey (LV)
					
Truffle (TF)	Cement (BT)	Acer (AH)	Corsica-Pear (KB)	Acacia (AK)	Norway Maple (SA)
					
Walnut (NB)	Beech (BU)	Oak-Basalt (EB)			

*Please note:
The pictured colour palette may
slightly deviate from colour and structure
of the original material!*

SW500



SW600



SW700



Laundry Trolleys SW

The laundry trolleys presented in the following section serve to streamline your work processes and offer the highest degree of hygienic safety.

All laundry trolleys consist of

- a buffed stainless steel frame with continuously welded pipe joints
- a cabinet made of special sheets and worktops, quality class E1 with double-sided melamine resin coating and a 2 mm thick ABS-safety edge foil
- 4 castor wheels Ø 125mm, with 2 brakes

In every series, the laundry trolleys are also available with a device to connect our linen collectors. They are equipped with an adhesive ring which is suitable for textile linen bags, linen nets and laundry bags.

Series SW600, SW700, SW900, SW1000 and SW2000 are also available with one or two drawers or a compartment that can be closed by roller blinds.

Linking Feature

Laundry trolleys with a linking feature to connect linen collectors are becoming increasingly popular. In two easy steps the linen collector can be attached or removed from the trolley. Our linking feature is designed such that linen collectors can be added freely to laundry trolleys already in use.

Standard colour: light grey (LG)

All available colours can be found in the colour palette on page 10.

! Appropriate accessories such as disinfection dispensers, glove boxes or writing plates can be found from page 71.

Article Number	Total Dimensions BxDxH mm
SW500	635 x 540 x 1080
SW600	635 x 590 x 1080
SW700	635 x 590 x 1275
SW800	490 x 580 x 1080
SW900	490 x 580 x 1080
SW1000	635 x 590 x 1275
SW2000	635 x 660 x 1010

Overview Laundry Trolleys

SW800



SW900



SW1000




SW2000

Image with optional railing



Linking Feature K

available for all laundry trolleys



Side compartments SF50

available for SW600/SW700



DEM

Lid with magnet
available for all laundry trolleys




Laundry Trolley with Integrated Collector System

SW620-KS70 **6** : model series
 2 : number of adhesive rings
 0 : number of drawers in the laundry trolley
 KS70 : waste collector at the front


front view with
2 drawers




2 adhesive rings




3 adhesive rings




4 adhesive rings




2 adhesive rings and KS70



3 adhesive rings and KS70



4 adhesive rings and KS70




Laundry Trolley with Coupled Collector System

SW620K-KS70 **6** : model series
 2 : number of adhesive rings
 0 : number of drawers in the laundry trolley
 K : linking Feature for collector system
 KS70 : waste collector at the front


front view with
2 drawers




2 adhesive rings




3 adhesive rings




4 adhesive rings




2 adhesive rings and KS70



3 adhesive rings and KS70



4 adhesive rings and KS70



Laundry Trolleys SW500



This series has open storage shelves for clean linen that can be accessed from three sides of the trolley. Standard colour is light grey.

The SW500 can be combined with all of our collectors (optional 2, 3 or 4 adhesive rings) - SW520, SW530 or SW540.

There are also systems with adapted linking features available - SW520K, SW530K, SW540K.

Additionally, the waste collector „Genius KS70“ (p. 82) can be mounted at the front of the laundry trolleys.

SW500



SW500K



SW500-W5



SW520



SW530

Image with optional accessories
KDS and TVD



Article Number	Total Dimensions BxDxH mm
SW500	635 x 540 x 1080
SW520	600 x 950 x 1080
SW530	600 x 1330 x 1080
SW540	600 x 1530 x 1080
SW500K	635 x 555 x 1080

Laundry trolleys of the series SW600 have a closed cabinet for clean linen (dimensions inside: 515 x 495 x 650 mm, BxDxH). The two storage shelves are adjustable in height and can be taken out completely. Doors that fully swing out allow easy access to the stored clean linen.

Also here, linking features can attach any collector either in a fixed or removable fashion.

All models can be provided with the waste collector „Genius KS70“ (page 82).

SW600



SW600K



SW620K

Image with optional accessories KDS and TVD



SW600-KS70



SW620



SW620-KS70



SW630



Article Number	Total Dimensions BxDxH mm
SW600	635 x 590 x 1080
SW620	605 x 1030 x 1080
SW630	605 x 1370 x 1080
SW640	605 x 1565 x 1080
SW600K	635 x 605 x 1080

Laundry Trolleys SW700

Series SW700 laundry trolleys have a closed cabinet for clean linen (dimension inside, B x D x H: 515 x 495 x 870 mm). The two storage shelves are adjustable in height and can be taken out completely. Doors fully swing out allow easy access to the stored clean linen.

Additional linking features can attach any collector either in a fixed or removable way. Furthermore, laundry trolleys of this series can also be equipped with either one or two drawers. All models can be provided with the waste collector „Genius KS70“ (page 82).

SW700



SW701K

Image with optional drawer with ZY



SW702K



SW720



SW721

Image with optional accessories
KDS
and
SPA2



SW732

Image with optional accessories
KDS and TVD



Article Number	Total Dimensions BxDxH mm
SW700	635 x 590 x 1275
SW701	635 x 590 x 1275
SW702	635 x 590 x 1275
SW720	605 x 1030 x 1275
SW730	605 x 1370 x 1275
SW740	605 x 1565 x 1275
SW700K	635 x 605 x 1275

Series SW800 models are distinguished by their reduced width of only 490 mm. Here the adhesive rings are located in a row, one behind the other. As with models of the SW500 series, this series has open storage shelves for clean linen that can be accessed from three sides of the trolley.

Linking features can attach almost any collector either in a fixed or removable fashion.

Models of this series can be provided with the waste collector „Genius KS70“ (page 82).

SW800



SW800K



SW810K-KS70

Image with optional accessory KDS



SW810

Image with optional accessory KDS



SW820

Image with optional accessories KDS and TVD



Article Number	Total Dimensions BxDxH mm
SW800	490 x 580 x 1080
SW810	490 x 980 x 1080
SW820	490 x 1300 x 1080
SW830	490 x 1595 x 1080
SW840	490 x 1890 x 1080
SW800K	490 x 595 x 1080

Laundry Trolleys SW900

SW900 are also distinguished by their reduced width of only 490 mm. The adhesive rings are located in a row, one behind the other. The difference from the SW800 series is the addition of a closed cabinet for clean linen and an access door.

The two storage shelves are adjustable in height and can be

taken out completely. Linking features can attach almost any collector either in a fixed or removable fashion.

Models of this series can be provided with the waste collector „Genius KS70“ (p. 82).

SW900



SW900K



SW910-KS70

Image with optional accessories
KDS and TVD



SW920K

Image with optional accessories
KDS and TVD



SW930

Image with optional accessories
KDS and
TVD



Article Number	Total Dimensions BxDxH mm
SW900	490 x 580 x 1080
SW910	490 x 1020 x 1080
SW920	490 x 1230 x 1080
SW930	490 x 1560 x 1080
SW940	490 x 1825 x 1080
SW900K	490 x 590 x 1080

Laundry Trolleys SW1000

Laundry trolleys of this series are built similar to the SW700 series. However, instead of a drawer system trolleys have a utensil rack, which can be closed and locked by a roller blind.

Linking features can attach almost any collector either in a fixed or removable fashion. Models of this series can be provided with the waste collector „Genius KS70“ (p. 82).

SW1000



SW1000K



SW1020

Image with optional accessory KDS



SW1010K

Image with optional accessories KDS, TVD and ZY



SW1020K-KS70

Image with optional accessories DMS5, TS5, KDS and TVD



Article Number	Total Dimensions BxDxH mm
SW1000	635 x 590 x 1275
SW1020	605 x 1030 x 1275
SW1030	605 x 1370 x 1275
SW1040	605 x 1565 x 1275
SW1000K	605 x 605 x 1275

Laundry Trolleys SW2000

Different from our other models and series, the SW2000 has a complete front that can be closed and locked by roller blind. The storage shelves are adjustable in height and can be taken out completely.

Linking features can attach almost any collector either in a fixed or removable fashion. Models of this series can be provided with the waste collector „Genius KS70“ (p. 82).

SW2000



SW2000K

Image with optional railing



SW2020



SW2020

Image with optional accessories KDS and TVD



SW2020K

Image with optional accessories KDS and TVD



Article Number	Total Dimensions BxDxH mm
SW2000	635 x 660 x 1010
SW2020	600 x 1075 x 1010
SW2030	600 x 1385 x 1010
SW2040	600 x 1505 x 1010
SW2000K	635 x 690 x 1010

List of Abbreviations

ABS	ISO-module, diff. vers.	FM	shelf, diff. vers.
ABSK	ISO basket	FP	base plate
ABSTOP	ABS synthetic top	FS8100	tea cart
ADEF	tray for defibrillator	GAE	stainless steel closeable dustbin, diff. vers.
AHS	fold-out protection for S70/S70L	GB	stainless steel grid shelf, various sizes
AK	stainless steel pull-out clothing rail	GF	spring-base trolley
ASD/ASK	stopping device for wire-netted and plastic baskets	GK	wire-netted basket, diff. vers.
ASM	stopping device for ABS- and PC-modules	GK2 - GK4	rubber toggle for linen nets, diff. vers.
ASPC4305-N	stopping device for PC-module PC4305-N	GSE	waste collector, diff. vers.
ATE	working and storage table, diff. vers.	GM	closed module trolley, diff. vers.
BS	push bar, diff. vers.	GS	aluminium cabinet trolley
BSML	B-Smart-Lock, electronic locking system	GTC/GTW/GTWK	transport trolley, diff. vers.
BSR	push bar for rack trolley	GTR2	skeleton container trolley
BTW01	trolley for abdominal bandages	H	textile protective cover for sorting racks, diff. vers.
CETOS	Phoenixcare-functional trolley, diff. vers.	help	fire extinguisher system
CI-...	OP equipment (shelves, cabinet, stools, tables)	HMF	magnetic holder mounted on chassis
CW	computer trolley, diff. vers.	HSB	glove box mounting, diff. vers.
D	set wire loops for drill machine	HW01	hook trolley
DD	drill machine	HWW	linen bags, warm- or hot water soluble
DDK	drill machine with chain	HZS	radiant heater for baby changing trolley
DEM	lid with magnet for all laundry trolleys	HZS-H	stainless steel mounting for radiant heaters
DESI	disinfection station	IAB	surgical instrument container with stainless steel fitting
DESICENTER	disinfection centre	IE	mobile instruments table, diff. vers.
DFST	door holder central, for hinged-folding door	IF	drip stand / IV-pole, diff. vers.
DFT	hinged-folding door	IFH	drip stand holder / IV-pole holder, diff. vers.
DIBE	white toggle with ring, diff. vers.	IFS	accessories for IV pole
DMS	disinfectant dispenser, diff. vers.	IPW	iPad-trolleys
DMSV	disinfectant closure plate, diff. vers.	IPH001	iPad mounting wall mounting
EAK	tubular frame, diff. vers.	K	linking feature
EBHES	stainless steel container for euro-rack	KB	synthetic waste bin, diff. vers.
EC/EE/EL	pedal dust bin, various sizes	KBW	synthetic base-tub, diff. vers.
EDVW	IT trolley, diff. vers.	KDH/KDS	loose plastic lid, diff. vers.
EDW	ergonomic documentation trolley	KF	clamping springs
EK	standardised label holders for rack systems	KI	lab coat stand
ENDO61000	endoscopic trolley	KKE, KKR	catheter tube, diff. vers.
ENDOHR	endoscopic bracket	KKL	brackets for border poles
EPD	stainless steel pump with dosing limitation for disinfectant dispenser	KNEBEL-ROT	toggle
ES	euro-rack, diff. vers.	KSR	combined sorting trolleys
EWK	stainless steel recyclables cabinet	KW01/KW02	catheter trolley, diff. vers.
F4/F6	compartment divider, various sizes	KW4/KW5	combined waste collector
FHL/FHR	F-mounting, mounted left/right	KW90	plastic tube, diff. vers.
FLEXISTORE	storage system FLEXISTORE	LAW	laptop-trolley, diff. vers.
		LGP	bag cassette for waste collector
		LTS	divider lengthwise, diff. vers.

List of abbreviations



MB01	mangle rack	SB4	bedpans holder
MB	Medibox	SCHUETTE	bins, diff. vers.
MBH	Medibox mounting euro-rack model	SEALSAFE	waste sealing unit
MBHS	Medibox mounting (set)	SEG/SEM/SEV	multi-purpose trolley
MEDI	medicine trolley, diff. vers.	SF50	Side compartments
MH/MH-VESA	monitor holder, diff. vers.	SFH01	oxygen-cylinder holder
MOD	module trolley series, diff. vers.	SH	"S"-hooks
MT1, MT2	furniture safe, diff. vers.	SLAZ450	heavy load extension drawer
MRSA	MRSA-trolley, diff. vers.	SLIDE	hygienic waste bin „SLIDE“, diff. vers.
MTW	girder support	SM	drawer module, diff. vers.
MUBH	mounting for mouth protection box, diff. vers.	SPA	extendable working surface, diff. vers.
N	linen net, diff. vers.	SR	sorting rack trolley, diff. vers.
NA	stainless steel waste collector "Compact", diff. vers.	SR400	adhesive ring, diff. vers.
NKV	file and ward round trolley, diff. vers.	SRW	shoe rack, diff. vers.
NOVOVAC	vacuum device NOVOVAC	SSR	shock protection disc
NRS01	kidney dish	SST01	bowl retainer
NSW	system trolley modula, diff. vers.	ST	standing racks base unit
NW	stainless steel linen collector "Compact" and „Compact Plus“, diff. vers.	STE	front wall trolley
NWE	recyclables collector Series 2000	SW	laundry trolley, diff. vers.
NWK	stainless steel recyclables cabinet Series 3000	TBW/TW	tray trolley, diff. vers.
NWKA	aluminium recyclables cabinet	TE	wall mounting, diff. vers.
OMW	open module trolley, diff. vers.	TEXT-MRSA	Sticker MRSA
OPT	operation room step	TH	T-mounting
P	corner posts, diff. vers.	TK	transport box, diff. vers.
PACTOSAFE	waste sealing unit, Pactosafe	TPS01	transponder lock
PAW	patient care and nursing trolley, diff. vers.	TRE	stainless steel transport trolley, diff. vers.
PC	ISO-module, diff. vers.	TRE-AP	hinged shelf for TRE
PCT	divider for ISO-module, diff. vers.	TS	drip pan, diff. vers.
PHA4ES	brochure holder for euro-rack	TSE	swab dispenser
PHE	stainless steel coat hooks	TVD	pedal with dumping function, diff. vers.
PHK	plastic coat hooks	TW1/TW2	linen bag, diff. vers.
PSS	power supply system	UF	urine bottles wall mounting, diff. vers.
QTS	slotted rail holders, diff. vers.	UTBES	implement basket
R	rack trolley, diff. vers.	WKK	waste collector with colored synthetic container
RG	border grid, diff. vers.	WW1/WW2	stainless steel undercarriage with plastic tube
RGE	tubular frame, diff. vers.	WW	baby changing trolley, diff. vers.
RIFIX	directional locking	WICKELAUFLAGE	baby's changing pad, diff. vers.
RLS	slotted rail holder	WSR	shoe rack for wall mounting
RP01	reanimation panel	WSU	stainless steel wall shelf
RSB, RSD	border poles, diff. vers.	WTW	linen bag, diff. vers.
RUWI	roll-up utility trolley	ZE	central lock
RW01	room cleaning trolley	ZFS	master adjusting device
S40/S70	waste collector „Genius“	ZS	combination lock
		ZSE	room recyclable Collector, diff. vers.
		ZY	cylindrical lock



Our goal is to provide you products and solutions that meet the latest design, technology and functional requirements.
Therefore, we have different processes to improve our products.

As part of these processes, there are optimizations as well as modifications or changes in production
that are relevant to your products.

Product changes are a sign of improvement and advancing manufacturing, but may cause variations
in color, style and texture of the products.

Errors and omissions excepted.



MERANTO OÜ

Siduri 3
11313 Tallinn/Estonia Tel.:
+372 50 181 04
sales@merianto.com
info@merianto.com
www.merianto.com
www.meriantomedical.com