


M-N Inner Cabinet for the Secure Storage of Psychotropic and Controlled Medicines

The Merianto Medical N-locker is designed for the separate and secure storage of psychotropic medicines, controlled substances, and other pharmaceuticals requiring enhanced security. The N-locker is installed inside a larger lockable medicine cabinet. This double-locking arrangement provides increased security and helps prevent unauthorized access to medicines. For enhanced security, we recommend using different locking systems for the outer medicine cabinet and the N-locker. The N-cabinet is fixed to the medicine cabinet structure using a minimum of six fixing points.


M-N1 locker

Dim: L270x350, H270mm
 1 adjustable glass shelf




M-N2 locker

Dim: L550x340, H390mm
 2 adjustable glass shelves




M-N3 locker

Dim: L410x210, H390mm
 2 adjustable glass shelves



M-N4 locker

Dim: L375x280, H375mm
 2 adjustable glass shelves



Product Features :

- welded steel construction, finished in white RAL 9003
- door panel located entirely within the cabinet body for improved resistance against forced entry
- height-adjustable internal glass shelves
- various locking options available (iLOQ, Abloy, RFID)
- N-lockers are dimensioned to fit a wide range of Merianto Medical medicine cabinets