


## MEDICINE CABINET MLK-65

Merianto's medicine cabinet **MLK-65** has been developed based on the ISO modular system taking into consideration the medical workers' need for a universal medicine cabinet. The lockable medicine cabinet has two sections and contains a separate metal safe with its own lock. The cabinet's interior can be fitted according to the user's need using ISO modular system 600 x 400 mm baskets or glass shelves.



The medicine cabinet is made from quality melamine. The facade can be a classical door (left/right) or a space conserving roller door. In the middle of the cabinet there is a pull-out auxiliary surface. The load capacity of shelves, drawers and baskets is 20 kg.



Standard options include five different finishings:



The standard solution for the medicine cabinet MLK-65 consists of:

- 19 mm thick melamine shell and doors
- Accessory hook at the cabinet top
- Clear Plexiglas pockets with height adjusting dividers on the upper door inner surface
- Pull-out auxiliary surface on ball-bearing runners
- Four glass shelves with adjustable height in the upper part
- A white RAL9010 lockable safe in the upper part
- ISO 600x400mm baskets (set choices from SET1-SET4) with dividers and label holders in the lower part
- LED lighting in the cabinet's upper part
- All doors have locks. The locks can be placed in series so that they are compatible with cabinets within the same room, floor or other desired zone. The base choice is a mechanical furniture lock, but electronic locks with codes or smart cards are additional options.
- Silicone magnets on doors – soft closing
- Strong metal foot frame with adjustable ends, white RAL9010 colored.

## MEDICINE CABINET MLK-65

### Lower cabinet basket selection (basket depth mm)



SET 1

50  
+  
100  
+  
100  
+  
100  
+  
100



SET 4

100  
+  
100  
+  
100  
+  
200



SET 2

50  
+  
100  
+  
200  
+  
200



SET 5

100  
+  
100  
+  
100  
+  
100  
+  
100



SET 3

50  
+  
50  
+  
50  
+  
50  
+  
100  
+  
100



SET 6

200  
+  
200  
+  
200